

# National Human Exposure Assessment Survey (NHEXAS)

## *Region 5 Study*

## Quality Systems and Implementation Plan for Human Exposure Assessment

Research Triangle Institute  
Research Triangle Park, NC 27079

Cooperative Agreement CR 821902

**Standard Operating Procedure**

**NHX/SOP-171-006**

**Title:** Maintenance of Perkin Elmer (PE) Model ZL5100 PC Atomic  
Absorption Spectrometer: Graphite Furnace

**Source:** Research Triangle Institute

U.S. Environmental Protection Agency  
Office of Research and Development  
Human Exposure & Atmospheric Sciences Division  
Human Exposure Research Branch

**Notice:** The U.S. Environmental Protection Agency (EPA), through its Office of Research and Development (ORD), partially funded and collaborated in the research described here. This protocol is part of the Quality Systems Implementation Plan (QSIP) that was reviewed by the EPA and approved for use in this demonstration/scoping study. Mention of trade names or commercial products does not constitute endorsement or recommendation by EPA for use.

**TITLE:** STANDARD OPERATING PROCEDURE FOR MAINTENANCE OF PERKIN  
ELMER (PE) MODEL ZL5100 PC ATOMIC ABSORPTION SPECTROMETER:  
GRAPHITE FURNACE

**SOURCE:** Research Triangle Institute  
Post Office Box 12194  
Analytical and Chemical Sciences  
Research Triangle Park, NC 27709-2194

**AUTHOR(s):**

Michelle Lang Date: 6/1/95  
\_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_ Date: \_\_\_\_\_  
\_\_\_\_\_ Date: \_\_\_\_\_

**APPROVED BY:**

Principal Investigator: E. Pellegrini Date: 6/2/95  
QA Officer: D.J. Smith Date: 6/7/95

<b><u>STATUS:</u></b>	IN PROGRESS:	<input type="checkbox"/>
	DRAFT:	<input type="checkbox"/>
	FINAL VERSION:	<input checked="" type="checkbox"/>

**REVISIONS:**

No.	Date	No.	Date
0	†	6	
1		7	
2		8	
3		9	
4		10	
5		11	

† Effective date of this version is the date of the last approval signature;  
revision 0 is the original version.

MAINTENANCE OF PERKIN ELMER (PE) MODEL ZL5100 PC  
ATOMIC ABSORPTION SPECTROMETER: GRAPHITE FURNACE

TABLE OF CONTENTS

<u>Section</u>		<u>Page</u>
<a href="#">1.0</a>	<a href="#">Introduction</a> .....	3
<a href="#">2.0</a>	<a href="#">Routine Maintenance Items</a> .....	3
<a href="#">3.0</a>	<a href="#">Non-Routine Maintenance Items</a> .....	3
<a href="#">4.0</a>	<a href="#">Documentation</a> .....	4
<a href="#">5.0</a>	<a href="#">Corrective Action</a> .....	4
<a href="#">6.0</a>	<a href="#">Sources/References</a> .....	4

## 1.0 INTRODUCTION

Maintenance practices for the PE Model ZL5100 PC Atomic Absorption Spectrometer graphite furnace operation are described below. They should be performed daily by the instrument operator during periods of analysis and recorded in the instrument log book. The PE User's Manual is kept near the instrument for reference.

## 2.0 ROUTINE MAINTENANCE ITEMS

### 2.1 Atomizer Cell

- Inspect graphite components for charring, pitting or other signs of degradation. Replace if necessary.
- Inspect quartz windows and o-rings; clean with methanol, ethanol, isopropanol, or other suitable solvent if film or particles present.

### 2.2 Controlled Temperature Furnace

- Check for adequate water level in recirculating chiller, and check for adequate flow.
- Check argon pressure and flow. Replace gas cylinder as needed.

### 2.3 A/S 70 Autosampler

- Check for adequate clean water (or appropriate rinsing solution, e.g., Triton-x, 0.5% HNO<sub>3</sub>) in rinsing reservoir.
- Empty waste bottle if necessary.
- Check integrity of injection tip (it should be clean, unclogged, unbent) and re-align sampler arm if necessary for accurate sample injection.

## 3.0 NON-ROUTINE MAINTENANCE ITEMS

### 3.1 Atomizer Cell

- The lab manager should remove and examine graphite contact rings annually. They should be cleaned with appropriate solvent and replaced if they are excessively charred or corroded.

- The lab manager should inspect the fume extraction system, including the graphite tip, tubing and motor, annually or in the event of system failure.

#### 4.0 DOCUMENTATION

The Instrument Logbook contains the daily record of use and conditions. Completion of routine maintenance is noted in this log book daily. Non-routine maintenance and repairs are recorded in a separate maintenance log section, with entries describing symptoms and repairs performed, and date of service.

#### 5.0 CORRECTIVE ACTION

If routine cleaning or replacement of consumable supplies described in Section 2.0 does not render the instrument operational, the lab manager should be consulted for further corrective actions. A Perkin Elmer service representative or RTI mechanical or electronics staff will be consulted as necessary.

#### 6.0 SOURCES/REFERENCES

1. Cantle, John Edward, ed. Atomic Absorption Spectrometry, New York, Elsevier, 1982.
2. Perkin Elmer Corp., Model ZL5100 PC Atomic Absorption Spectrometer Users Manual, Norwalk, CT, 1987.
3. Perkin Elmer Corp., Model ZL5100 PC Software, Version 6, Supplemental Information, Norwalk, CT, 1991.
4. Rains, Theodore C. "Atomic absorption spectrophotometry" in Water Analysis, Vol. II, Part 2, R. A. Minear and L. H. Keith, eds., Academic Press, 1984, pp. 61-110.
5. Willard, H. H.; Merritt, L. L.; Dean, J. A. Instrumental Methods of Analysis, New York, D. van Nostrand Co., fifth ed., 1974.